## ACTA DE INSPECCION

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fixed to the motor shaft, rotates within the indexer channel encoder. Signals from the encoder provide the microprocessor with information as to which channel is selected. The transfer tube/adaptor optopair (11), comprising an infrared light source and a detector sensitive to light, gives information about whether a transfer tube/adaptor is connected to the selected channel. The locking ring optopair (13) detects if the transfer tubes/adaptors are locked in the indexer channels.

The microprocessor subsequently locates the check cable tip and the source in the required positions within the applicators by counting steps made by the stepper motor, starting from a reference point, during the out-drives. The reference optopair (12), located at that point, is used to determine when the check cable tip or the tip of the source capsule is passing it. For the out-drives of the cables this moment is reached as soon as the light path between the optopair is obstructed. For both cables the number of steps required to restore the light path (pass the optopair) during the in-drive is compared with the number of steps made during the out-drive. Within certain limits, these numbers must correspond.

The degree of friction between the cables and the applicator is dynamically measured during the out and in-drives by determining the number of steps that the shaft lags behind the number of stepper motor pulses given by the microprocessor.

The source or check cable is retracted if an unacceptable difference is detected during the out-drive. If this happens during the in- or out-drive of the source, the system activates the emergency stop circuit if the normal retraction using the stepper motor fails.

**The emergency stop motor (4) is coupled to the shaft of the source stepper motor via an electrically powered clutch. When the source is not in the safe and the emergency stop circuit is activated it removes electrical power from the source and check cable stepper motors, while turning on the emergency stop motor. The latter then retracts the source until the wire-in switch (3) is actuated (indicating that the source is in the safe) or a maximum time period has elapsed.**

The emergency stop motor is powered by a DC power supply with battery backup allowing it to operate independently from the mains supply. The maximum pull it exerts exceeds the pull of the stepper motor by a factor of five. This force is still well below the level that would damage the cable or the weld between the cable and the source capsule.

## 6.2.2 Description of the Indexer

The indexer stepper motor positions the channel selection tube of the indexer (see Fig. 6-24) at the programmed channel. Attached to the tube is the indexer disc, which has a different combination of holes for each channel.

For the 18-channel indexer, the disc has two test positions (19 and 20). For the 30-channel indexer, the disc also has two test positions (between 1 and 2, and between 2 and 3). The holes are sensed with five opto-pairs. During the self test and before a fraction/pulse, the system checks for correct operation of the indexer.

The signals from the opto-pairs are routed to the opto-pair interface board, where they are conditioned for transfer to the power supply and control board.

Channel	M4/	M3/	M2/	M1/	M0/
1	1	1	1	1	0
2	1	1	1	0	1
3	1	1	1	0	0
4	1	1	0	1	1
5	1	1	0	1	0
6	1	1	0	0	1
7	1	1	0	0	0
8	1	0	1	1	1
9	1	0	1	1	0
10	1	0	1	0	1
11	1	0	1	0	0

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